

### Technical Data - Multiwave 7000

Technical Specifications					
	Rack 24	Rack 18	Rack 9	Rack 6	Rack 5
Number of vials	24	18	9	6	5
Quartz vial	Yes	Yes	Yes	Yes	Yes
PTFE-TFM vial (HF-resistant)	Yes	Yes	Yes	No	Yes
Glass vial (disposable)	Yes	Yes	Yes	Yes	Yes
Rec. filling volume	5 mL	10 mL	25 mL	40 mL	55 mL
Sample amounts (organic)	up to 0.2 g	up to 1 g	up to 2 g	up to 3 g	up to 4 g
Sealed quartz vessel	No	Yes	No	Yes	No
Applications	All kinds of samples can be digested with all racks. Choose your rack according to whether you want to run with pressure-sealed vials or sealed quartz vessels and according to the required sample amount, sample volume, or reaction mixture.				